




Please Read Carefully Before Applying ProCollar to your pet



(1) To Secure ProCollar in place:

- Make sure that your pet's *everyday collar* goes through the inner ring straps of the ProCollar.
- Securely fasten pet's everyday collar.

(2) Inflate ProCollar via air valve (*pinch & blow*).

(3) Connect Velcro strap on the ProCollar.

- ✓ *For larger breeds with longer limbs, or breeds in between sizes, always size up!*
- ✓ *Always monitor a pet post-surgery to ensure collar is working properly. Some pets may be able to maneuver around any collar.*